



SpecCheck Secondary Gel Standard Kit - Fluoride, 0-2.0 mg/L F

Proizvod #: 2712500

EUR Cijena: Kontaktirajte nas
Otprema unutar 2 tjedna

Solid, ready for use Gel-Standards for instrument verification in the lab and in the field.

Set contains four vials filled with gels that simulate the test colour at various concentrations.

Specifikacije

Koncentracija:	0 – 2,00 mg/L F
Naziv metode:	SPADNS
Uređaj:	All except DR/3000, DR 900 and DR/800 Series