



Quality control standard hardness, low range, 500 mL

Proizvod #:	2833449
EUR Cijena:	Kontaktirajte nas
Dostupno	

Ready to use - no dilutions required. Mixed standard includes: Total Hardness: 100 mg/L as CaCO₃. Calcium Hardness: 50 mg/L as CaCO₃.

Specifikacije

Koncentracija:	Calcium Hardness 50 mg/L CaCO ₃
Opis:	Low Range Hardness, Mixed Parameter QC Standard
Platforma:	Otopina