



Quality control standard oxygen demand , pk/16 - 10 mL Voluette ampule

Proizvod #:	2833510
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Ready to use - no dilutions required. Oxygen Demand Quality Control Standard. BOD - 396 mg/L. COD - 617 mg/L. TOC - 245 mg/L.

Specifikacije

Koncentracija:	BOD – 396 mg/L
Opis:	Oxygen Demand, Mixed Parameter QC Standard
Platforma:	Otopina