



Quality control standard metals for drinking water, high range, 500 mL

Proizvod #: 2833649

EUR Cijena: Kontaktirajte nas

Nije prikazan datum isporuke

Ready to use - no dilutions required. Mixed standard includes: Copper: 2.5 mg/L Cu, Iron: 1.5 mg/L Fe, Manganese: 5 mg/L Mn. 500 mL Bottle

Specifikacije

Koncentracija:	Mangan: 5,00 mg/L Mn
Opis:	Metals Drinking Water (High Range) Mixed Parameter QC Standard
Platforma:	Otopina