



Quality control standard metals for drinking water, low range, 500 mL

Proizvod #: 2833749
EUR Cijena: Kontaktirajte nas
Dostupno

Ready to use - no dilutions required. Mixed standard includes: Copper: 1 mg/L Cu, Iron: 0.3 mg/L Fe, Manganese: 0.1 mg/L mg/lMn. 500 mL Bottle.

Specifikacije

Koncentracija:	Mangan: 0,1 mg/L Mn
Opis:	Metals Drinking Water (Low Range) Mixed Parameter QC Standard
Platforma:	Otopina