



Sodium Standard Solution, 10 mg/L Na, 1 L

Proizvod #:	2835153
EUR Cijena:	Kontaktirajte nas
Dostupno	

NaCl Standard, 10 mg/L Na, 1 L. For Polymetron 9240 & 9245 analyser.

Specifikacije

Čuvanje uređaja:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Koncentracija:	10 mg/L
Parametar:	Sodium
Platforma:	Otopina
Uređaj:	9245 & 9240
Volumen/ veličina paketa:	1 L 1 L