



Bottle, BOD, 300 mL, 24/pk, serialized: 1-24

Proizvod #: 2898700
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included Serialized bottles have silk-screened, sequential numbers (1-24) and include a bar code.

Specifikacije

Body Material:	Glass
Kapacitet:	300 mL
Materijal:	Borosilicate Glass
Package size UOM:	24/pk, Serialized: 1 - 24
Volumen/ veličina paketa:	24 /pk