



Be Right™



Bottle, storage, autoclavable, wide mouth, 125 mL pk/12

Proizvod #:	2939103
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Polyethylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121°C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Specifikacije

Body Material:	Plastic
Kapacitet:	125 mL
Materijal:	PP
Package size UOM:	12pk
Volument/ veličina paketa:	12 /pk