



Fixed Volume, Micropipette, 20 μ L

Proizvod #:	2946003
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and color-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Specifikacije

Kapacitet:	20 μ L
Točnost:	0,5 %
Volumen/ veličina paketa:	1 each