



Be Right™



pH paper, 1.0 - 11.0 pH range

Proizvod #:	39133
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Provides simple, fast, economical method of pH measurement. Sharp color change for each 0.5 pH unit. Convenient for sample pH adjustments. Recommended for sewage water, not finished water. Five 4.6-meter (15') rolls per package.

Specifikacije

Mjerni rang:	1.0 - 11.0 pH units
Volument/ veličina paketa:	5 4.6-meter (15') rolls per package