



## Syringe, Luer-Lok Tip, individually packaged, 3 cc

|                              |                   |
|------------------------------|-------------------|
| <b>Proizvod #:</b>           | <b>4321300</b>    |
| EUR Cijena:                  | Kontaktirajte nas |
| Nije prikazan datum isporuke |                   |

Disposable syringes feature a Luer-Lok tip suited for connection with syringe filters or filter holders. Marked graduations. 3 cc, individually packaged, 1/pk.

---

### Specifikacije

|                           |        |
|---------------------------|--------|
| Kapacitet:                | 3 mL   |
| Volumen/ veličina paketa: | 1 each |