



Nitrate Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

Proizvod #:

4456369

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Otprema unutar 2 tjedna



Ionic Strength Adjustor (ISA) for Nitrate determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Nitrate level readings. Eliminate most interferences. Pack of 100 powder pillows.

Specifikacije

Broj testova:	100
Koncentracija:	ISA
Mjerni rang:	0,05 - 14.000 mg/L NO ₃ ⁻ -N
Naziv metode:	Ion Selective Electrode
Opis:	Nitrate Ionic Strength Adjuster Powder Pillows
Package settings:	Powder Pillows
Parametar:	Nitrate
Volumen/ veličina paketa:	100 /pk