



Beaker, Griffin, Low Form, Glass, 600 mL

Proizvod #:	50052
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifikacije

Jedinica:	1 each
Kapacitet:	600 mL
Materijal:	Glass, Griffin, Low Form
Materijal tijela mjerača:	Glass