



Beaker, griffin, low form, glass, 2000 mL

Proizvod #:	50054
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Specifikacije

Jedinica:	1 each
Kapacitet:	2000 mL
Materijal:	Glass, Griffin, Low Form
Materijal tijela mjerača:	Glass