



Beaker, griffin, low form, glass, 100 mL, 12/pk

Proizvod #:	50072H
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Specifikacije

Jedinica:	12 /pk
Kapacitet:	100 mL
Materijal:	Glass, Griffin, Low Form
Materijal tijela mjerača:	Glass