



Be Right™



Flask, filtering, 500 mL

Proizvod #:	54649
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Glass filtering flask with side-arm at 90° angle. Heavy-walled with graduated capacity scale. Designed to ASTM E1406, Type II, Class I requirements. Hose connection accepts 5/16" I.D. flexible tubing.

Specifikacije

Jedinica:	1 each
Kapacitet:	500 mL
Materijal:	Glass
Materijal tijela mjerača:	Glass