



Portable HQD rugged field case for one or two rugged probes with extended cable length

Proizvod #:

8505600

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Rugged, Convenient Storage for Field Use

Rugged Field Case features custom cut foam inserts to hold portable meter, a maximum of 2 rugged probes with up to 10 meters of cable, and accessories in place.

All-in-One durable rugged field cases provide convenience and include space for sample collection containers, meter and rugged probes

Durable, die-cut foam inserts hold meter and probes securely in place

Specifikacije

Dimenzije (V x Š x D):

420 mm x 458 mm x 216 mm