



Recharge kit of 4 pre-filled cartridges with premounted 2956A membranes for A1100 oxygen sensors with 33051-xT cap only

Proizvod #:

2956A-AT05

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Otprema unutar 2 tjedna



Includes o-rings and cleaning tools. Special low concentration electrolyte.

Specifikacije

Čuvanje uređaja:

Store in a dry place at room temperature

Parametar:

O₂

Uređaj:

Orbisphere A1100-S00T O₂ sensor, Orbisphere GA2x00-S00T O₂ sensor