



## Replacement Sample Cup for SL1000

|                    |                   |
|--------------------|-------------------|
| <b>Proizvod #:</b> | <b>9418100</b>    |
| EUR Cijena:        | Kontaktirajte nas |
| Dostupno           |                   |

Replacement Sample Cup for SL1000

---

### Specifikacije

|                           |                            |
|---------------------------|----------------------------|
| Što je u pretincu:        | One Replacement Sample Cup |
| Volumen/ veličina paketa: | 1 each                     |

---

### Što je u pretincu

One Replacement Sample Cup