



## Formaldehyde cuvette test - ISO 12460, 0.5-10 mg/L H<sub>2</sub>CO, 25 tests

**Proizvod #:**

**LCK425**

EUR Cijena:

Kontaktirajte nas

Dostupno

### **LCK cuvette tests: Minimum amounts of reagents, maximum safety.**

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

**Precise and reliable measurement values**

**Maximum safety for users, thanks to the closed system and low amounts of reagents**

**Convenient and error-free dosing of the reagents**

**Barcode label for automatic recognition in the photometer**

**Differentiation between tests and measuring ranges by means of colour-coding**

---

### **Specifikacije**

Broj testova:	25
Čuvanje uređaja:	15 - 25 °C
Mjerni rang:	0,5 - 10 mg/L H <sub>2</sub> CO
Naziv metode:	ISO12460
Parametar:	Formaldehyde
Platforma:	LCK
Uređaj:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000