



Be Right™



pH pufer otopina, pH 10,01, 500 mL

Proizvod #:	2283649.99
EUR Cijena:	Kontaktirajte nas
Dostupno	

Pufer otopina za kvalitetu, označena plavom bojom, NIST sljediva.

Specifikacije

Kodirani po boji?:	Da
Koncentracija:	10.01 pH $\pm$ 0.02 pH при 25° C
Package settings:	Boca
Parametar:	pH
pH vrijednost-i:	10,01 pH $\pm$ 0,02 pH pri 25°C
Platforma:	Otopina
Volument/ veličina paketa:	500 mL