



Be Right™



Otopina za punjenje, referentna, 1 M KNO_3 , 59 mL

Proizvod #: 28564026
HRK Cijena: Kontaktirajte nas
Dostupno

Otopina za punjenje referentnog dijela elektrode, kalijev nitrat (KNO_3), 1 M, boca od 59 mL

Specifikacije

Formula:	KNO_3 1 M
Opis:	MTC306
Solution Type:	Ponovno punjenje
Volument/ veličina paketa:	59 mL

Što je u pretincu

Otopina, 59 mL