



pH Chemkey Reagents (box of 25)

Proizvod #:	9759000
EUR Cijena:	Kontaktirajte nas
Dostupno	

Use with the SL1000 Portable Parallel Analyser (PPA) for dramatically streamlined water quality testing.

pH Chemkeys, range 6.3 - 9.0. For SL1000 Portable Parallel Analyser (PPA).

Avoid the hassle of pH probes and calibration in the field with our easy to use colorimetric pH Chemkey. The Hach SL1000 Portable Parallel Analyser (PPA) performs the same tests with less than half the manual steps. Get highly accurate results, with less opportunity for errors, in a fraction of the time. Up to six parameters, tested in parallel. The SL1000 is a portable colorimeter that provides faster testing by using Chemkey reagents instead of powder pillows. All chemicals and processes are entirely contained inside the Chemkey.

Current SL1000 instrument customers who want to test pH need to perform a software update. That software can be located on the Downloads tab of the SL1000 instrument page.

Specifikacije

Mjerni rang:	pH 6.3 - 9.0
Parametar:	pH
Platforma:	Chemkeys
Uređaj:	SL1000
Volument/ veličina paketa:	Kutija od 25 komada