



## Acid solution for EZ4004 Total Alkalinity Analyser, 0.1 N HCl (50-2500 mg/L Range), 2 L

**Proizvod #:**

APPC4004-03-EU **Opasno po zdravlje**

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke



### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Specifikacije

Koncentracija:	0.1 N HCl
Parametar:	Total Alkalinity
Volument/ veličina paketa:	2 L