



## Buffer solution for EZ1016 Total Hardness Analyser

**Proizvod #:**

**APPC1016-02-EU** **Opasno po zdravlje**

**EUR Cijena:**

**Kontaktirajte nas**

**Nije prikazan datum isporuke**



## Ready-to-use Reagents for EZ Series Online Analysers

### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

## Specifikacije

Parametar:

Total Hardness