



Cylinder, Graduated, Mixing, Glass, 50 mL $\pm$ 0.4 mL, 1.0 mL divisions, glass stopper #16

Proizvod #: 104592
EUR Cijena: Kontaktirajte nas

Tall form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Specifikacije

Certificiranja:	Class B
Kapacitet:	50 mL
Materijal tijela mjerača:	Glass
Mjerna metoda:	To Contain
Volumen/ veličina paketa:	1 each
Vrsta:	Jednometrički, miješanje