



**Be Right™**



## **Orthophosphate Standard Solution, 15 mg/L PO<sub>4</sub>-P, 500 mL**

<b>Proizvod #:</b>	<b>LCW1056</b>
EUR Cijena:	Kontaktirajte nas
Dostupno	

Standard Solution for use with the NP6000sc Phosphate Analyser.

---

### **Specifikacije**

Koncentracija:	15 mg/L PO <sub>4</sub> -P
Parametar:	Orthophosphate
Uređaj:	NP6000sc
Volument/ veličina paketa:	500 mL