



Be Right™



General Boiler Water Test Kit, Model AL-94 MG-L (Drop Count Titration)

Proizvod #:

222403-EU

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Used for controlling boiler and cooling water treatment and checking corrosion, this test kit is an economical solution for the small boiler operator who needs to monitor alkalinity, chloride, hardness, iron, and pH.

Uses simple, economical drop count titration and color disc methods

Kit contains all required reagents and apparatus to perform analysis in the field

Uses drop count titration and color disc methods

Specifikacije

| | |
|------------|---|
| Model: | AL-94 MG-L |
| Parametar: | Alkalinity, chloride, hardness, iron, and pH (4-10) |
| Platforma: | Disk sa bojama-Color Disc |
| Primjena: | Boiler Water |