



Hardness (Calcium and Total) Reagent Set, Digital Titrator

Proizvod #:

2272100-EU

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



For Calcium and Total Hardness determination by EDTA Digital Titrator titration. Ranges: 10-160 and 100-4000 mg/L CaCO₃. Reagent Set contains: Buffer Solution (Hardness 1), CalVer 2 Calcium Indicator Powder Pillows, ManVer 2 Hardness Indicator Powder Pillows, Potassium Hydroxide Solution (8 N), EDTA Titration Cartridges (1 each) (0.800 and 0.0800 M), and 5.25N sulfuric acid.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Mjerni rang:	10 - 160 mg/L & 100-4000 mg/L CaCO ₃
Naziv metode:	EDTA Titration
Parametar:	Hardness (Calcium and Total)
Posebna značajka:	Reagent set, Digital Titrator

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)