



Be Right™



NitraVer 5 Nitrate Reagent Powder Pillows, 5 mL, pk/100

Proizvod #:

1403599-EU

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



Powder reagent in sealed packages for the determination of Nitrogen-Nitrate up to 50 mg/l NO₃-N.

For high range nitrate determination by the NitraVer 5 cadmium reduction method. Range: up to 50 mg/L as nitrate-N. Sample Size: 5 mL. Pack of 100 powder pillows.

Specifikacije

Mjerni rang:	5 mL
Naziv metode:	Cadmium Reduction
Platforma:	Praškasti reagensi PP
Volument/ veličina paketa:	100 /pk

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)