



Ammonia Ionic Strength Adjustor (ISA), Solution, 500 mL

Proizvod #:

2824349-EU

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke



Ionic Strength Adjustor (ISA) for Ammonia determinations by Ion Selective Electrode (ISE) method. Eliminate most interferences. Solution, 500 mL bottle.

Specifikacije

Koncentracija:	ISA
Package settings:	Bottle
Parametar:	Ammonium
Volumen/ veličina paketa:	500 mL