



Be Right™



Test kit, nitrate-nitrite, model NI-12

Proizvod #:
EUR Cijena:
Nije prikazan datum isporuke

1408100-EU
Kontaktirajte nas

Opasno po zdravlje



Simple, low-cost color disc method is the most accurate visual test kit method.

Specifikacije

Broj testova:	100
Mjerni rang:	0 - 40 mg/L NO_3^- -N
Mjerni rang (2):	0,01 - 0,40 mg/L NO_2^- -N
Model:	NI-12
Naziv metode:	Color disc/Cadmium reduction, Diazotization
Parametar:	Nitrogen, Nitrate/Nitrite - NO_3^- -N/ NO_2^- -N
Parameter:	Nitrogen, Nitrate/Nitrite
Platforma:	Disk sa bojama-Color Disc
Smallest increment:	1 0.01
Standard range:	0 - 40 mg/L NO_3^- -N; 0.01 - 0.40 mg/L NO_2^- -N
Test Kit Model:	NI-12**
Vrsta kućišta:	D